

Addition to the mounting guide.

Thank you for purchasing the 194 Bench Wall.
As it is not released officially, there are still some work in progress.

1. With this version of the product you will receive 2× 4mm washers.

You will have to mount the bench through the mounting interface shown by the circles here (pink and blue). The 4mm washers should be added to the backside between the wall and the bench frame (blue circle). This is to make sure the mounted bench doesn't have a few degrees forward tilt and stay level.

2. A little force is required when mounting the seat to the steel frame. Mount one end first with all 4× M10 screws (green). Then one should bend/force the tube a bit while another inserts the M10 screws. (Yellow)

NOTICE: That the pink and blue mounting interfaces doesn't allow for the same diameter of threaded rods to go through.

Check mounting guide for step by step method.

